

Urologicals



Urinary Drainage Bag - Sterile

Reorder No.	Quantity
4271	20/Cs

One (1) each 2000 ml latex free drainage bag equipped with:

- Anti Reflux Valve
- Drip Chamber
- Air Vent
- Universal Hanger
- Urine Sampling Adapter

Sterility of the interior of drainage bag and fluid pathway is guaranteed if the connector cap is in place, drain tube is closed and unit is undamaged.

Competitor Amsino - Ref # AS32200

Advantage Urinary Drainage Bag - Sterile

Reorder No.	Quantity
4270	20/Cs

One (1) each 2000 ml latex free drainage bag equipped with:

- Anti Reflux Valve
- Air Vent
- Universal Hanger
- Urine Sampling Adapter

Sterility of the interior of drainage bag and fluid pathway is guaranteed if the connector cap is in place, drain tube is closed and unit is undamaged.

Competitor Amsino - Ref # UB2000TV

