



Suction Canister



Rest Easy with Dynarex Respiratory Products

- Durable plastic container system designed for collecting and retaining aspirated materials
- Float shut-off with built-in filter increases flow, eliminating damaging backflow to the suction regulator
- Foam guard and float valve ultrasonically welded in place, preventing premature shutoff
- Port caps attached to lid, which saves time
- Patient port designed with 20° angle (1200cc version) to prevent kinked tubing – no separate adapter required
- Universal design easily fits wall mounts



1200cc or 800cc Canister



Reorder No.	Description	Quantity
4670	1200cc Suction Canister (Hi-Flow) w/Lid	40/cs
4671	1200cc Suction Canister Kit (Hi-Flow) w/Lid with 3/16" dia. Tubing – 18" and 6'	30/cs
4672	800cc Suction Canister (Hi-Flow) w/Lid	24/cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/ sizes currently not shown.



Also see our website for Suction Tubing (#4680-4689)

© 2012 Dynarex Corporation

For more information contact your Dynarex sales representative 888-396-2739

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarexassets.com