



Yankauer Suction Handle

Breathe Easy with Respiratory Care Products from Dynarex

- **Sterile suction tool for removing particulates and viscous material**
- **Bulb tip to minimize tissue trauma**
- **Ribbed connector for secure attachment to suction tubing**
- **Blue-tinted, transparent material with slip-resistant handle**
- **Smooth inner surface for maximum flow**
- **Vented and non-vented models**
- **Latex free**

Reorder No.	Description	Quantity
4690	Yankauer Suction Handle – vented	50/cs
4691	Yankauer Suction Handle – non-vented	50/cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.



Also see our website for our Suction Canisters (#4670-4672)

© 2012 Dynarex Corporation

For more information contact your Dynarex sales representative
888-396-2739

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarexassets.com