



Dressing Change Tray – Sterile



Convenience in a Complete and Sterile Kit

Dynarex Dressing Change Tray is a convenient sterile kit for aseptic dressing changes. All necessary items are packed in a single, disposable tray, eliminating the need to maintain inventory on multiple items and assembling them for a procedure.

- **Sterile**
- **Latex free**
- **Components arranged in procedural sequence for convenience**
- **Cost-effective and disposable**

Contents:

- | | |
|--|---|
| (1) Instructions for use | (1) Waterproof drape, 45 x 34cm |
| (1) Wound care tray | (1) Dressing change label |
| (1) Non-woven face mask | (1) PVC gloves, size M, powder free, white, with wallet |
| (1) Non-woven split sponge, 2" x 2", 4 ply | (1) Tape measure, 36" |
| (1) Non-adherent sponge, 2" x 2" | (3) PVP swabsticks |
| (1) Gauze sponge, 2" x 2", 8 ply | (3) Alcohol swabsticks |
| (1) Gauze sponge, 4" x 4", 8 ply | (1) 3M™ tape, 3/4" x 18" |
| | (1) 3M™ Tegaderm transparent film dressing, 4" x 4-3/4" |



Room No. _____
Name _____
Dr. _____



Reorder No.	Description	Quantity
4705	Dressing Change Tray – Sterile	20/cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

Visit our website to see our Minor Laceration Tray (#4523) and our Staple/Suture Removal Kits (#4521-4522)



© 2012 Dynarex Corporation

For more information contact your Dynarex sales representative
888-396-2739

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarexassets.com