




Tillotson True Advantage Nitrile Exam Gloves



Discover the True Advantage

Manufactured using a patented non-sulfur based curing process which connects the polymer molecules unlike the curing methods of other nitrile gloves. Allows the polymer more freedom of movement, allowing the glove to conform to the hand rather than constrict it.

- The first incredibly soft nitrile glove that is chemo approved and improves upon the feel of latex
- Non-Latex Powder Free Exam Glove with polyglide lining



Reorder No.	Description	Quantity
6501	Nitrile Exam Glove - XS	10/100/Cs
6502	Nitrile Exam Glove - S	10/100/Cs
6503	Nitrile Exam Glove - M	10/100/Cs
6504	Nitrile Exam Glove - L	10/100/Cs
6505	Nitrile Exam Glove - XL	10/100/Cs

Chemotherapy Approved

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

See our website for our Pure Advantage Sterile Surgical Gloves (9160-9185)



Already mobile? Scan your smart phone here for more information.



Find us on:



For more information contact your Dynarex sales representative 888-396-2739