



SensiLance™ Pressure & Button Activated Safety Lancets

For Painless Sampling

SensiLance™ Pressure Activated Lancets

- For capillary blood samples
- Ready to use - no loading needed
- Needle automatically retracts
- Latex Free

SensiLance™ Button Activated Lancets

- Tri-bevel tip for virtually painless sampling
- Needle automatically retracts
- For capillary blood samples
- No lancing device required
- Latex Free



Already mobile?
Scan your
smart phone
here for more
information.



Reorder No.	Description	Quantity
7112	21 gauge - 2.2 mm Depth, Pressure Activated	10/100/cs
7114	26 gauge - 1.8 mm Depth, Pressure Activated	10/100/cs
7115	28 gauge - 1.8 mm Depth, Pressure Activated	10/100/cs
7122	21 gauge - 2.4 mm Depth, Button Activated	10/100/cs
7124	26 gauge - 1.8 mm Depth, Button Activated	10/100/cs
7125	28 gauge - 1.8 mm Depth, Button Activated	10/100/cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.
7112-7125 HCPC: A4259



See our website for
our complete line of
exam gloves

For more information
contact your Dynarex sales representative
888-396-2739

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2012 Dynarex Corporation