

About EZ-ACCESS

EZ-ACCESS is the world's leading manufacturer of high quality aluminum ramp systems. We are a family-owned and operated company with over 30 years of experience in designing and manufacturing access solutions. Our products are made in the USA and are known for their safety, durability, and workmanship.

In addition to the versatile GATEWAY™ ramp, we manufacture the Passport® Vertical Platform Lift and a variety of portable ramps and other mobility accessories for both commercial and residential use.

EZ-ACCESS can make any type of home, regardless of age or style, accessible. Our products provide flexible short- and long-term solutions for a variety of needs.

Through an extensive network of local authorized dealers, EZ-ACCESS serves the entire United States with professional products and quick installation.

For additional information, please contact us today or visit ezaccess.com.



customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com



GATEWAY™

Solid Surface Portable Ramp

Life changes.

Be ready for the unexpected with help from EZ-ACCESS.



GATEWAY Portable Wheelchair Ramp

Designed for safety and engineered for reliability, the GATEWAY, a solid surface portable wheelchair ramp, is a smart alternative to costly wooden ramps. Available with or without handrails.

Visit ezaccess.com for product specifications

//// GATEWAY™ Solid Surface Portable Ramp

Designed for safety and engineered for reliability, the GATEWAY is a smart alternative to costly wooden ramps. When mobility needs extend beyond the scope of a conventional portable ramp, the durable GATEWAY ramp is your answer.

The GATEWAY Advantage

- The safe and durable home access solution.
- Solid, seam-free surface. The 1-piece design means no assembly* or movable parts; just set the ramp in the desired position, and it's ready to use!
- Made of lightweight, yet heavy-duty aluminum.
- Designed for the rigors of repeated, everyday use in all weather conditions. Can be easily moved or relocated as needs change.
- Makes a great semi-permanent ramp for indoor/outdoor use.
- Made of durable aluminum for years of maintenance-free use. Will not rust, rot, or deteriorate from the outdoor elements.
- Permanent (extruded) slip-resistant surface allows superb traction.
- Available with or without handrails.
 - Handrails are welded for extra strength and durability.
 - Feature a brushed, architectural finish.
- Designed for scooters, wheelchairs, and walkers.
- Three-year limited warranty.
- Made in the USA.

* Handrails require minimal assembly.



The EZ-ACCESS® Choice

Safety is our top priority, and we are committed to making the process of choosing the best ramp solution as worry-free as possible. With EZ-ACCESS, you can be confident that your accessibility needs will be resolved safely and quickly.

The Aluminum Advantage

The decision on whether to install an aluminum accessibility solution or build a wood ramp can be a difficult one. However, aluminum ramps have many advantages over their wooden counterparts.

- Virtually maintenance free
- Will not rot or rust
- Easy assembly with minimal tools
- Environmentally friendly
- An asset that can be reused, recycled, or moved to another location if needed

Quick Specs

- Welded aluminum construction
- Ramps feature a standard width of 36 inches
- Permanent (extruded) slip-resistant surface
- Weight capacity of 850 pounds
- Available in 3-, 4-, 5-, 6-, 7-, 8-, 9-, & 10-foot lengths

//// Call today for more information ////

customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com