

## 2-Point Slings

- One piece sling for use with 2-point spreader.
- Provides support from shoulders to knees.
- Used for transferring, weighing or toileting.
- For use with chain set, Item GF133-S-C (not included with sling, sold separately).
- Canvas slings are green. Nylon slings are white.
- Designed to meet the requirements of HCPCS Code: E0621.
- 6 month limited warranty.

### 220 lb. Weight Capacity

#### Canvas Fabric

GF112-C-LC	Without Commode Opening	1 ea
GF113-C-LC	With Commode Opening	1 ea


#### Nylon Fabric

GF112-N-LC	Without Commode Opening	1 ea
GF113-N-LC	With Commode Opening	1 ea

### 400 lb. Weight Capacity

#### Heavy Duty Canvas Fabric

GF114-C-LC	Without Commode Opening	1 ea
GF115-C-LC	With Commode Opening	1 ea



Sling Type	Body Width	Body Length	Commode Opening Width	Commode Opening Length
Without Commode Opening	26"	41 ½"	NA	NA
With Commode Opening	26"	41 ½"	6 ¾"	10 ⅝"



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**www.grahamfield.com**

**Fax: 1-800-726-0601**

GF0700193RevA07