

Digital Lift Scale

- Digital scale easily attaches to Lumex Patient Lifts.
- Completely self-contained.
- Combines the latest microprocessor and load cell technology to provide accurate and repeatable weight data with a calibrated accuracy of 0.1%.
- Measures pounds or kilograms.
- Maximum weight capacity 600 lb.
- Comes complete with attachment hardware.
- 1 year limited warranty.

LF1009A Fits LF1030, LF1030CH and LF1040 1 ea

DSC200 Fits LF1050 and LF1090 1 ea

The same digital scale as described above comes complete with spreader, arm, and hardware for attachment to the LF1075 SURE-LIFT P3 Patient Lift .

LF1007 1 ea



LF1009A



DSC200



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com