



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Narcotics Cabinet, Large, Stainless Steel Line, Single Door/Double Lock #2732

Specifications



Construction:

- Heavy-Duty, 20 gauge stainless steel

Features:

- True double locking protection with 2 high security pick-resistant tubular locks
- Locks are keyed separately, with four keys per lock included
- Unique lock sequences— all cabinets are keyed differently
- One fixed shelf and one adjustable and removable shelf with the choice of 4 locations
- Full length stainless steel pinned door hinges prevent pin removal
- Predrilled holes in the back for mounting on a wall, cabinet or other solid surface (wall mounting hardware included)
- Narcotics cabinet is delivered ready to use (door handle must be customer installed)
- Twelve year warranty against manufacturer defect

Dimensions:

- 29.5" H x 23.5" W x 10.5"D (74.9 cm H x 59.7 cm W x 27.7 cm D)
- Net weight: 35 pounds (15.9kg)
- Approx. shipping weight: 40 pounds (18 kg)



One lock and key set includes black identifier