

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini Line Phlebotomy Cart Package with Key Lock #3245KPB



Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (1) 6-inch (15.2 cm)
- (1) 9-inch (22.9 cm)

Sixteen available color choices: *Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Features

- Wrap-around Stabilizing Bumper
- 5-inch (12.7 cm) premium plate casters
- Cabinet has 24 vertical inches (60.96 cm) drawer space
- · Key-locking system with 2 keys
- · Pre-threaded accessory mounting holes
- · Ball-bearing full extension drawer slides
- Removable plastic top
- · Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish available in 16 colors

Accessories Included:

- 2-tier raised back rail system (MINI-RBR)
- Extended Push Handle
- Single Glove Box Holder (AL2362)
- Aluminum Waste Container with Lid (AL680436)
- Two 6-bin tilt bin organizer (683805)

Other Mini Line Phlebotomy Cart Options:

- Locking Bemis ® Sharps Holder with 5 quart Container (684801)
- Upgrade to 3 Box Glove Holder (680422UPG)
- Chart Holder (680417)

Specifications

Construction:

· Painted steel cabinet and drawers

Internal Drawer Dimensions:

• 15" W x 16.75" D (38.1cm x 42.5cm)

Assembled Cart Dimensions

- 52.3" H x 27.5" W x 21.5 "D (133 cm H x 70 cm W x 54.6 cm D)
- Height to Work Surface 33.8"
- Cart net weight: 108 lbs (49.1 kg)
- Approximate shipping weight: 133 pounds (28.2 kg)