



650 Ford Street
Colorado Springs, CO 80915
Toll: (800) 433-4064 * Ph: (719) 637-0300
Fax: (719) 597-8273 * www.harloff.com

Classic Line, Six Drawer, Emergency Cart, Standard Package #6400

Specifications



Construction:

- Sturdy, uni-body steel construction

Features:

- Stainless steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Stainless steel top rail / push handle
- Full wrap-around bumper on separate stabilizing frame
- Five-inch (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered fully assembled
- Twelve year warranty

Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions

- 44.5" H x 32" W x 22" D
(113 cm H x 81.3 W x 55.9 D)
- Cart weight: 131 lbs
- Approximate shipping weight: 151 pounds