



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Nine Drawer, Pediatric Resuscitation Cart, Super-Stat Specialty Package #6411PEC



Assembled Cart Dimensions

- 44.5" H x 32" W x 22" D
(113 cm H x 81.3 W x 55.9 D)
- Cart weight: 165 lbs (74.8 kg)
- Approx shipping weight: 200 pounds (91 kg)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Drawer Configuration:

- (8) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 50 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Wrap-around bumper on separate stabilizing frame
 - The 6411PEC features a super extended width bumper that is designed to enclose and protect all included accessories such as the suction unit shelf, etc
- Includes a pack of multi-colored labels for the end user to apply to identify drawers
- Wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty

Accessories Included:

- Cardiac board and brackets (680406)
- Three utility hooks (680402)
- IV pole (680401)
- Hospital grade 6-socket outlet strip with 15' cord and cord wrap (680404)
- Stainless steel defibrillator shelf 20" x 12" (680321)
- Suction unit shelf (680308)
- Oxygen tank holder - D/E size cylinders (680403)
- Chart holder (680416)
- Locking sharps container, 5qt sharps, and glove box holder (684811)
- Two large utility baskets (681501)
- One small utility basket (680431)