



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Six Drawer Workstation (Tall) with Key Lock, Standard Package #6450



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Two aluminum side rails / push handles
- Full wrap-around vinyl padded bumper on separate steel frame
- Twelve year warranty
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (4) – 3"
- (1) – 6"
- (1) – 12"

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D
(105.4 cm H x 81.3 W x 55.9 D)
- Cart net weight: 140 lbs
- Approximate shipping weight: 165 lbs