



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 6 Drawer Workstation (Tall) with Key Lock, Deluxe Package #6456



### Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

### Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D  
(105.4 cm H x 81.3 W x 55.9 D)
- With Accessories: 55.9" H x 44" W x 22.1" D
- Cart net weight: 171 lbs
- Approximate shipping weight: 195 lbs

### Sixteen available color choices:

#### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

## Specifications

### Construction:

- Sturdy, uni-body steel construction

### Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Full wrap-around bumper on separate steel frame
- Durable, hard-baked powder coat finish
- Twelve year warranty

### Specialty Package Accessories Included:

- Raised back rail system with two horizontal rails (68163)
- Four aluminum side rails / push handles (68161)
- 6-Compartment tilt bin organizer (683805)
- Tape and label dispenser (680415)
- Catheter box, 18", with rail clip (680413)
- Large utility basket (680407)
- 6-Socket hospital grade outlet with 15-foot cord and cord wrap (680404)
- Plastic waste container, 3 gallon, with bracket and rail clip (684801, 683801)
- Large sharps container and bracket with rail clip (684802, 683801)
- 3-inch exchange trays (2) (681501)
- 6-inch exchange trays (2) (681502)
- 12-inch Adjustable drawer divider set (680503)

### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)