



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Difficult Airway Cart, 6 Drawer Workstation (Tall) with Key Lock, Specialty Package #6457



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Two aluminum side rails / push handles
- Full wrap-around bumper on separate steel frame
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty

Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Assembled Cart Dimensions

- 58.4" H x 40" W x 25" D
(148.3 cm H x 101.6 W x 63.5 D)
With full height of adjustable shelf: 62.6" / 159 cm
- Cart net weight: 190lbs (86.1 kg)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Difficult Airway Accessory Package Includes:

- Extended width bumper (68030)
- Stainless steel adjustable overhead shelf with accessory bar (68032)
- Oxygen tank holder, D/E cylinders, stainless steel (680403)
- Hospital grade 6-socket outlet strip with 15' cord and cord wrap (680404)
- Scope cabinet – side mounted; steel with clear plexi-glass door, key locking. Holds four scopes (680350)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)