



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Eight Drawer Workstation, Key Lock, Standard Package #6610

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with two keys
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Aluminum side rail / push handle
- Full wrap-around bumper on separate stabilizing frame
- Twelve year warranty
- Latex free
- Durable, hard-baked powder coat finish



Drawer Configuration:

- (6) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D
(105.4 cm H x 81.3 W x 55.9 D)
- Cart net weight: 145 lbs
- Approximate shipping weight: 165 lbs