



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Five Drawer Aluminum Mini Line Anesthesia Cart with Extended Stabilizing Bumper and 5" Wheels, with Basic Electronic Pushbutton Lock #AL3245E

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints

Specifications

Construction:

- Painted aluminum cabinet and drawers

Features:

- Lightweight aluminum construction
- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 37" H x 19.5" W x 25" D
(94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 64 pounds (29 kg)
- Approximate shipping weight: 84 pounds (38.1 kg)



**Color shown is Sand with Hammertone Green Drawers
Additional Colors Available**