



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## **Five Drawer Aluminum Mini Line Anesthesia Cart with Extended Stabilizing Bumper and 5" Wheels, with Key Lock #AL3245K**

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints

### **Specifications**

#### **Construction:**

- Painted aluminum cabinet and drawers

#### **Features:**

- Lightweight aluminum construction
- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key lock with two keys
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

#### **Drawer Configuration:**

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### **Internal Drawer Dimensions:**

- 12.6" W x 16.75" D (32cm x 42.5cm)

#### **Assembled Cart Dimensions:**

- 37" H x 19.5" W x 25" D  
(94.0 cm H x 49.5 cm W x 63.5 cm D)
- Cart net weight: 64 pounds (29 kg)
- Approximate shipping weight: 84 pounds (38.1 kg)



#### **Sixteen available color choices:**

##### **\*Antimicrobial Finish Available**

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White