

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Aluminum Universal Line Medium 5 Drawer Procedure Cart, with Basic Electronic Pushbutton Lock #AL808E5

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.



Color Shown in Light Blue.
Additional Colors Available.

Specifications

Construction:

Aluminum construction

Features:

- Internal pullout shelf on right side
- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- 26.0" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- · Cart is delivered wheels-on and mostly assembled
- Pre-drilled threaded holes for standard accessories
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year warranty on manufactured parts

Drawer Configuration:

- (3) 3.25" (8.3 cm)
- (1) 6.5" (16.5 cm)
- (1) 9.75" (24.8 cm)

Internal Drawer Dimensions:

• 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 41" H x 31" W x 24.75" D
- (104.1 cm x 78.7 cm x 62.9 cm)
- Cart weight: 105 pounds (47.8 kg)
- Approximate shipping weight: 135 pounds (62 kg)