



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Universal Line Super 7 Drawer Cart with Eletronic Lock AL810EKC7

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.

Specifications

Construction:

- All Aluminum Construction

Features:

- Electronic-locking cart with keypad access and BEST® key lock override with 2 keys
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- Allows for supervisor assignments; dual credential access available
- May add, edit, or delete users and user groups
- Battery-operated with six D-cell batteries; life of one year (minimum) at 25 to 50 uses/day
- Lock has an LED low battery indicator
- All aluminum, light-weight construction
- 32.5" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper with side extension for Oxygen Tank holder
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled hole locations for standard accessories – many are pre-threaded for ease of installation
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year manufacturer's warranty

Drawer Configuration:

- (5) – 3.25" (8.3 cm)
- (1) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Assembled Cart Dimensions:

- 48" H x 31" W x 24.75" D (121.9 cm H x 78.7 cm W x 62.9 cm D)
- Cart weight: 110 pounds (50 kg)
- Approximate shipping weight: 135 pounds (61.2 kg)

