

Standard Isolation Station ISO6550



Specifications

650 Ford Street

Colorado Springs, CO 80915

Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Construction:

Painted steel cabinet, drawers and shelves

Features:

- Silver ion paint coating with antimicrobial properties
- Unique open well design for visual storage of isolation supplies
- Well lid can be latched open
- Folding side shelves (2) with quick-release latch
- Built-in glove box holder (holds up to three glove boxes) and mask box holder
- Area for hand sanitizer/soap dispenser
- · Adjustable coat and hat rack
- · Hard Baked paint finish
- · Patent-pending design
- Ball-bearing, full extension drawer slides
- 5" ball bearing premium rubber casters, full swivel, 1 with toe-locking brake and 1 directional
- · Cart is delivered fully assembled

Thirteen available color choices:

Standard color choice: Cream Yellow Call for other choices

can for carer cricicos

Assembled Cart Dimensions

Open: 60.5" H x 46.63" W x 27.2" D (153.67 cm x 118.5 cm x 69 cm) With coat and hat rack extended 66.36" H Closed: 50.31" H x 25.68" W x 27.2" D (128 cm x 65.3 cm x 69 cm) Cart net weight: 102 pounds (46.4 kg)

Approx. shipping weight: 133 pounds (60.5 kg)

Open Well Dimensions:

• 13" H x 21.5" W x 22" D (33 cm H x 54.61 cm W x 55.88 cm D)

Internal Drawer Dimensions:

• 20"W x 17.25"D (50.8 cm W x 43.81 cm D)

Drawer Configuration:

- 1-6" Storage drawer
- 1-9" Storage drawer