

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# Classic Line Medication Cart for Multi-Dose Packaging with Key Lock #MD50-6K



#### Thirteen available color choices:

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray and Hammer Tone Green

## **Specifications**

#### **Construction:**

· Painted steel cabinet and drawers

### Features:

- Tall cabinet with key lock system
- Vinyl-padded full wrap-around steel bumper
- Five-inch sealed ball-bearing premium casters, all full swivel one with brake, one directional
- Internal locking narcotics box in 1/3 of bottom drawer
- Drop shelf
- · Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- · Certified "Latex Free"
- Five year warranty
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- · Hard-baked powder coat finish
- Available in Specialty Package. (MD50-6KSP)
  "SP" package includes sharps and waste container with brackets.

## **Drawer Configuration:**

- 5 4" Drawers with plastic drawer dividers (10 spaces per drawer)
- 1 10" Drawer with 1 adjustable divider and locking narcotics box

#### **Assembled Cart Dimensions:**

- 41.5" H x 32" W x 22" D (105 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 155 pounds (70.5 kg)
- Shipping weight: 180 pounds (81.8 kg)