



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Preconfigured Double Column Medical Storage Cart with Tambour Door MS8140-LAP

The MS8140-LAP is a preconfigured double column medical storage cart, with tambour door, designed for quick access to medical supplies, featuring the most innovative modular tray system on the market. The MS8140-LAP is designed to store supplies for surgical procedures related to the abdomen or pelvis through small incisions with the aid of a camera to inspect and diagnose a condition or to perform surgery.

The MS8140-LAP is preconfigured as follows: left column, top to bottom, 2ea 1/2" trays, 1ea 2" tray, 2ea 4" trays, 1ea 8" tray, 2ea 4" wired baskets, 3ea 7" wired baskets; right column, top to bottom, 2ea 1/2" trays, 1ea 2" tray, 2ea 4" trays, 1ea 8" tray, 2ea 4" wired baskets, 3ea 7" wired baskets.



Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains five hooks.

Standard Features:

- One half inch modular tray: Total of four trays that are 23.50 x 15.75 x 1/2 inch
- Two inch modular tray: Total of two trays that are 23.50 x 15.75 x 2 inches
- Four inch modular tray: Total of four trays that are 23.50 x 15.75 x 4 inches
- Four inch wired baskets: Total of four baskets that are 23.50 x 15.75 x 4 inches
- Seven inch wired baskets: Total of six baskets that are 23.50 x 15.75 x 7 inches
- Eight inch modular trays: Total of two trays that are 23.50 x 15.75 x 8 inches
- Tray dividers to organize width and depth of modular trays are sold separately
- Plastic side panels form a modular system consisting of 94 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 42" W x 29.75" D
(204.47 cm H x 106.68 cm W x 75.565 cm D)
- Cart net weight: 247 pounds (112.04 kg)
- Approximate shipping weight: 295 pounds (133.81 kg)