



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* [www.harloff.com](http://www.harloff.com)

## Double Wide Full Shell Catheter Procedure Cart with Tambour Door Designed for use in a Cath Lab #MS8140CATH

The MS8140CATH is a full shell, double wide Catheter Procedure Cart with tambour door, specifically designed for use in a cath lab to support the catheterization procedure



Close Up Of Cath. Slide Shelf

### Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

### Advantages and Benefits:

- The significant amount of catheter sizes and packaging techniques has made organizing and managing this critical area an extremely difficult task. Harloff offers a solution with a high density Catheter Storage System, which provides quick access and retrieval of catheters and related supplies.
- Equipped with adjustable Cath Slide Shelf, consisting of seven full extension catheter slides, each slide has five hooks, and each hook can hold up to 15 catheter packages.
- Each Cath Slide features a pair of ball bearing slides to withstand the pull force of removing catheters from hooks; full extension brings the catheters into full view for quick access and visual inventory of remaining stock.
- Cath Slide Shelf can be mounted in top position or 13 inches lower for storage above.
- The MS8140CATH stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access with our specially-engineered catheter.
- The maximum package length is 61 inches.

### Standard Features:

- Catheter Slide Shelf, featuring seven cath slides and five hooks per slide, can be mounted in a lower position to create storage above for over-sized items.
- Cath slides consist of a pair of ball bearing drawer slides to withstand the pull force of removing catheters from hooks.
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty

### Assembled Cart Dimensions:

- 80.5" H x 42" W x 27.9" D  
( 204.5 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 350 pounds (158 g)
- Approximate shipping weight: 480 pounds (217 kg)