



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Full Shell Medical Storage Cart Equipped with Five Shelves and Tambour Door #MS8140SH

The MS8140SH is a full shell medical storage cart equipped with five shelves and tambour door designed for quick access to medical instruments and supplies, featuring the most innovative modular tray system on the market. The MS8140SH can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.



Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- The MS8140SH has an open shell with five removable shelves. The open shell design offers 17% more storage capacity when compared with leading competitor's product.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- The shelves, with five dividers per shelf, are designed to improve visibility, accessibility and organization of sterile supply storage in the operating room.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Plastic side panels allow for additional shelves to be mounted.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains five hooks.

Standard Features:

- Five adjustable shelves with five adjustable dividers per shelf. Dividers can be mounted on ½ inch increments
- Plastic side panels form a modular system consisting of 47 insertion opportunities for mounting shelves
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Color: Sand, durable hard-baked hammer-beige powder coat textured finish
- Ships fully assembled with twelve year warranty

Assembled Cart Dimensions:

- 80.5" H x 42" W x 27.9" D
(204.5 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 275 pounds (124.7 kg)
- Approximate shipping weight: 405 pounds (183.7 kg)