



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Standard Line 360 Punch Card Medication Cart, Standard Package #SL360PC

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication cart with 360 punch card capacity, single width internal narcotics box
- BEST® brand lock on cabinet. Illinois lock on narcotics drawer
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Painted steel pull-out shelf
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL360PCSP) Included are:
#68640-P1 drawer tray with adjustable dividers, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea - 9.75" punch card drawer with dividers and locking narcotics box
- 2 ea – 9.75" punch card drawers with dividers

Internal Drawer Dimensions:

- 18" D x 20.5" W

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 26.5" W x 23.5" D
(116.2 cm x 67.3 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 137 lbs (62.2 kg)
- Approx shipping weight: 167 lbs (75.9 kg)