



650 Ford Street
Colorado Springs, CO 80915
Toll: (800)433-4064 * Ph: (719)637-0300
Fax: (719)597-8273 * www.harloff.com

Standard Line, Wide, Seven Drawer Anesthesia Cart, Key Lock, Standard Package #SL7KW



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Extra wide cart for your special storage needs
- Contemporary design with rounded corners
- BEST® lock key-locking system with 2 keys
- Pull-out shelf on right side of cart
- Full wrap-around bumper on stabilizing frame
- 5" ball-bearing casters; full swivel – one with brake, one tracking
- Ball-bearing, full extension drawer slides
- Removable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- Twelve year warranty
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (1) - 4" (10.2 cm)
- (4) - 3.25" (8.3 cm)
- (1) - 6.5" (16.5 cm)
- (1) - 9.75" (24.8 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 27.6" W x 17" D (70.1 cm x 43.2 cm)

Assembled Cart Dimensions

- 45.75" H x 34.2" W x 23.75" D
(116.2 cm H x 86.9 cm W x 70 cm D)
- Cart net weight: 137 pounds (62.2 kg)