



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## Wood Laminate 480 Capacity Punch Card Medication Cart, Standard Package #WL480PC-CM

This wood look mobile punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



### Specifications

#### Construction:

- Aluminum base and internal frame
- MDF panels covered with a special heat-pressed thermal foil laminate with a rich wood appearance

#### Features:


- Single column medication cart with 480 punch card capacity
- BEST® brand lock on cabinet
- Illinois lock on narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Painted steel pull-out shelf
- 4" ball-bearing premium rubber casters all full swivel
- Bottom skirting conceals casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

**Additional Wood Look Option Available,  
#WV540PC-CM Vinyl Coated Aluminum**



Cherry  
Mahogany  
(CM)

WL480PC-CM



Solara  
Oak  
(SO)

WL480PC-SO

### Color Options:

Color shown in Cherry Mahogany with a gray granite finish on top. Also available in Solara Oak #WL480PC-SO

Additional laminate color options available with minimum cabinet purchase and extended lead time.

### Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

### Assembled Cart Dimensions:

- 45" H x 35.5" W x 22.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 220 lbs (100 kg)
- Approx shipping weight: 290 lbs (131.8 kg)