

B. Model 1449-57
Shown with option:
3002-57 – Mat



C. 3002-57 Mat attached with VELCRO® brand fasteners



B. Series 1449
Dual Lift Powermatic® Mat
Platforms with Black Laminate Tops



- Will accommodate 2" thick reversible, removable mat, ordered separately. See #3002 below.
- Black laminate rimless subtop.
- Power height adjustment from 19" wheelchair access height to 28"H.
- Two direct drive synchronized lift mechanisms for unsurpassed power.
- "Hands-free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed

MODEL	WIDTH	LENGTH	HEIGHT	OPT.PAD
1449-46	4'	6'	19"-28"	3002-46
1449-47	4'	7'	19"-28"	3002-47
1449-57	5'	7'	19"-28"	3002-57
1449-68	6'	8'	19"-28"	3002-68

C. Series 3002
Platform or Floor Mats



- Nylon-reinforced vinyl mat with 2" thick 100 I.L.D. urethane foam padding.
- (4) handles for easy moving and hanging.

COLORS:

Reversible – Blue on one side, Gray on the other with a Gray 2" border.

MODEL	WIDTH	LENGTH	HEIGHT
3002-46	4'	6'	2"
3002-47	4'	7'	2"
3002-57	5'	7'	2"
3002-68	6'	8'	2"

D. Series 3042
Folding Floor Mats




- Accordion folding, 2' wide panels.
- 2" extra-firm foam.
- Nylon-reinforced vinyl cover.
- Easily joined to other mats with Velcro® brand fasteners on mat ends.
- Two-tone Blue and Gray.

MODEL	WIDTH	LENGTH	HEIGHT
3042-46	4'	6'	2"
3042-48	4'	8'	2"



C. Model 3002-47



D. Model 3042-46