



**A. Model 4431**



**B. Model 4440-LD**  
Shown with options:  
15 – Safety Rail  
091 – Pelvic Tilt

## A. Models 4431, 4436 Econo-Line™ Power Clinic Examination Table



### FEATURES:

- Dual direct-drive electrical lift actuators adjust height from 24" to 35".
- Convenient hand control operation.
- (2) door, (2) drawer storage cabinet.
- Includes 36" long (10) position adjustable backrest. (Must raise past 60° to lower.)
- 3 position adjustable 18" long legrest.
- Vinyl-clad stirrups adjust in length.
- Includes Paper Dispenser and Cutter.
- 110V, 60Hz, 1 Amp. Hospital Grade Plug.

**COLORS:** (FOLD OUT FRONT COVER)

**Laminate** – Folkstone Gray.

**Upholstery** – 17 standard vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4431	56"-74"*	27"	24"-35"
4436	Left Hand Table		

\*Fully extended

## B. Series 4440, 4445 Powermatic® Exam Tables

### FEATURES:

- Motorized lift from 26" to 38" high.
- "Hands-Free" foot control.
- Removable top pad with built-in pillow is attached with VELCRO® brand fasteners.
- Backrest is pneumatic Air-Spring with dual adjustment controls from 0 to 75°.
- Sliding 21" wide x 16" long leg rest.
- Concealed stirrups adjust in length and adjust laterally to (4) positions.
- (2) laminate drawers in front & side.
- Full width black plastic Drawer pulls.
- Black laminate base with levelers.
- Includes Paper Dispenser and Cutter.
- 110 V, 60Hz, 3.5 Amps.



**COLORS:** (FOLD OUT FRONT COVER)

**Laminate** – 927 Folkstone Gray or your choice of any high pressure laminate solid color or woodgrain.

**Upholstery** – 17 standard vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4440	55" to 71"	26"	26" - 38"
4440-LD (UL)	Same but less Duplex outlet		
4445	Left Hand Table		
4445-LD (UL)	Same but less Duplex outlet		

## OPTIONS FOR SERIES 4440, 4445

- \*04 – Safety Strap
- \*15 – (1) Safety Rail (Order [2] for both sides)
- \*17 – Utility Rails (both sides)
- \*089 – Motorized Power Backrest
- \*091 – Pelvic Tilt (Controls on both sides)
- \*097 – Stainless Steel Pan and Cutout

\*FACTORY INSTALLED

SEE PAGE 4