

A. Model 4643



see Stools
pages 30, 31

Common Features

TWO-IN-ONE EXAMINATION/TREATMENT TABLES — MODELS 4602, 4624, 4643:



- Hausmann Two-In-One Tables save space and costs by providing a readily adjustable examination and treatment table in a single, attractive and functional unit.
- Extended, the table is 72" long with a 14" adjustable headrest, a 17" adjustable footrest and a 41" center section.
- Adjustable length stirrups are included as standard equipment.
- 2" High-Density urethane foam top.
- Resist-All™ vinyl upholstery finish.

COLORS:

Laminate – Natural Oak, Wild Cherry, Folkstone Gray

Upholstery – 17 standard vinyl colors
(FOLD OUT FRONT COVER)

A. Model 4643

SPECIAL FEATURES:

- Fully enclosed all laminate cabinet.
- (2) hinged doors with 3mm edging.
- Concealed, self-closing Euro-hinges.
- Laminate interiors with (1) adjustable shelf.
- Black laminated plywood base.



MODEL	LENGTH	WIDTH	HEIGHT
4643	41" - 72"*	24"	31"

* Fully extended.

B. Model 4602

SPECIAL FEATURES:

- H-Brace is recessed into high pressure laminate table legs with unique 4-sided "lock-tite" joint.

MODEL	LENGTH	WIDTH	HEIGHT
4602	41" - 72"*	24"	31"

* Fully extended.

C. Model 4624

SPECIAL FEATURES:

- H-Brace with full width laminate shelf.
- High pressure laminate legs.

MODEL	LENGTH	WIDTH	HEIGHT
4624	41" - 72"*	24"	31"

* Fully extended

OPTIONS


FOR 4600 SERIES TABLES

- *04 – Table Safety Strap
- *19 – Slide Shelf
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- *99L – Left Hand Table Model

*FACTORY INSTALLED

SEE PAGE 4

B. Model 4602



C. Model 4624

