



A. Model 4746

Common Features

POWERMATIC® TREATMENT TABLES —
MODELS 4746, 4747:



- The Powermatic® electro-mechanical drive mechanism provides an infinite number of height positions.
- Patients can mount or dismount without a footstep. You can work faster and more efficiently without back strain or fatigue.
- 2" High-Density urethane foam vinyl top.
- "Hands Free" foot control operation.
- (2) laminate drawers on nylon rollers.
- Hinged door storage compartment.
- Includes Paper Dispenser and Cutter.
- Hospital grade plug.
- 110V, 60Hz, 3.5 Amps. U.L. Listed.
- New contoured black base plate with (4) levelers and satin aluminum trim.
- Shipped fully assembled.

COLORS:

Laminate – Folkstone Gray or choice of any high pressure laminate solid color or woodgrain.

Upholstery – 17 standard vinyl colors.

(FOLD OUT FRONT COVER)

A. Model 4746 Powermatic® Treatment Table with Gas-Spring Backrest

SPECIAL FEATURES:

- Gas-spring pneumatic 30" backrest elevates to 75° with lever control on one side.
- Power height adjusts from 26" to 38".

MODEL	LENGTH	WIDTH	HEIGHT
4746	76"	27"	26"-38"

Model 4747 (not shown) Powermatic® Treatment Table with Power Backrest

SPECIAL FEATURES:

- Same as Model #4746 except Power 30" backrest can be elevated automatically to 75° with foot control.

MODEL	LENGTH	WIDTH	HEIGHT
4747	76"	27"	26"-38"

Option: 15 – Safety Rail is mounted at mid-section to one side or both sides of table.



OPTIONS

FOR MODELS 4746, 4747

- *04 – Table Safety Strap
- *15 – (1) Safety Rail
(Order [2] #15 if both sides)
- 32 – Pillow, full size
- 35 – Pillow, small
- *99L – Left Hand Table

*FACTORY INSTALLED

SEE PAGE 4