



## SELECT FROM A VARIETY OF BARIATRIC PRODUCTS

4723



SEE PAGE 13

4416



SEE PAGE 19

4795



SEE PAGE 35

2038



SEE PAGE 34

6045



SEE PAGE 39

1429



SEE PAGE 43

6185



SEE PAGE 62

1386



SEE PAGE 71

1564



SEE PAGE 78

### C. Model 6053 3-Section Hi-Lo Multi-Position Treatment Table



#### FEATURES:

- Multi-position table can be used for a wide variety of therapy and treatment functions.
- 76"L x 27"W 3-section 1½" firm density gray vinyl top features contoured face/nose cutout in head section.
- Motorized height adjustment from 19" wheelchair height to 37"H.
- "Hands-free" Foot Control operation.
- Adjustable head section, 15" x 27", raises up to 45° and lowers down to 90° below horizontal.
- 22"L x 27"W middle seat section is stationary
- 39"L x 27"W Leg/Backrest section raises easily up to 73° with sturdy friction locking device and dual release handles.
- Top sections have welded steel support frames.
- Includes (4) 3" diameter swiveling casters with locking brakes.
- 110V, 60 Hz, 1.8 amps. U.L. Listed.

MODEL	LENGTH	WIDTH	HEIGHT
6053	76"	27"	19" - 37"

C. Model 6053

