

### A. MODEL A9066 TREATMENT TABLE WITH OPEN SHELF STORAGE



- Open shelf storage, accessible from both sides of the table.
- Center section has fixed shelves. Left and right sections have (1) adjustable shelf.
- Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- 10 Pro-Form™ vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9066	78"	30"	31"



A. MODEL A9066

### B. MODEL A9067 CABINET TREATMENT TABLE



- (6) Hinged doors on both sides of table.
- Storage accessible from both sides of table.
- Center section has fixed shelves. Left and right sections have (1) adjustable shelf.
- Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- 10 Pro-Form™ vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9067	78"	30"	31"



B. MODEL A9067

#### OPTIONS

MODELS A9066, A9067, A9073

\*69 - Gas Spring Backrest 0° to 75° (A9066, A9067 center fixed shelves move to head section.)

\*A946 - (1) Lock/Bolt for Double Doors

\*FACTORY INSTALLED



Model A9067 shown with doors open

### C. MODEL A9073 CABINET TABLE



- (3) Two-door cabinets with laminate interiors and adjustable shelf.
- Solid finished back.
- Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- 10 Pro-Form™ vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9073	78"	30"	31"



C. MODEL A9073