

A. Model A9098 Crank Hydraulic Taping Table

SPECIAL FEATURES:

- End open cabinet and supply shelf.

MODEL	LENGTH	WIDTH	HEIGHT
A9098	48"	30"	30" - 40"

Common Features:

CRANK TAPING TABLES – MODELS A9062, A9098:

- Crank hydraulic lift with fold-away handle allows height adjustment for back saving efficiency.
- Hydraulic cylinders are built into (4) black powder coated steel legs.
- 2" thick High-Density foam vinyl upholstered top.

COLORS:

Laminate – Natural Oak

Upholstery – 9 Pro-Form vinyl colors
(FOLD OUT FRONT COVER)



A. Model A9098



B. Model A9062



Shown as Treatment Table
with middle pad in place

B. Model A9062 Crank Hydraulic Taping/ Treatment Table

SPECIAL FEATURES:

- Two Person Taping Table converts into a full-length Treatment Table by adding middle pad.
- Fixed surface under middle pad is laminated for work use.
- Storage cabinet has (2) full extension Pull-Outs and (2) laminate drawers.
- 250 lbs. per seat for a combined capacity of 500 lbs.

MODEL	LENGTH	WIDTH	HEIGHT
A9062	78"	36"	30" - 40"

C. Models A9590, A9594 Pro-Lockers

- Dowel and glue laminated wood construction with 3mm high-impact edging.
- Foot Locker with vented storage and hinged laminate seat (upholstered on A9594) with hasp for padlock.
- Double prong clothes hooks.
- Chrome clothes rod on A9590 only.
- Quantity discount available for contract orders.

MODEL	WIDTH	DEPTH	HEIGHT
A9590	32"	22"	78"
A9594	18"	24"	78"

OPTIONS

MODELS 9590, 9594

- A9591 – Locking Security Door (for A9590)
- A9595 – Locking Security Door (for A9594)



C. Models A9590, A9594
Shown with optional locking
security door

