

Capacity: 330 lb / 150 kg



Advantages

- ✓ **6 -1/2" Dial Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **Low Cost**
This high quality scale provides accurate weight measurements at an affordable price
- ✓ **Weighs in LB or KG**
Dial displays shows weight in both pounds and kilograms providing versatility

Capacity: 330 lb / 150 kg

142KL Dial



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
142KL	2	20.13" x 12.5" x 8.5" / 1.24	14 lb	892076002572

Specifications

- Capacity: 330 lb / 150 kg
- Graduation: 1 lb / 0.5 kg
- 6-1/2" dial
- Platform Size: 9-7/16" (w) x 10-3/8" (d)
- 1 Year Limited Warranty