




Advantages

- Sanitary Plastic Weighing Surface**
 The easy-to-clean plastic tray is ideal for sterile environments
- Built-In Measuring Tape**
 Easy-to-read measuring tape allows user to measure both weight and height simultaneously
- Engraved Numerals**
 Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability
- Adjustable Unit Measurement Bar**
 The mechanical rotating poise bars allow for easy reading of weight in lb or kg

Health o meter[®]
P r o f e s s i o n a l



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
1522KL	1	20.75" x 25.25" x 11.5" / 3.49	35 lb	892076002176



0470

Specifications

- Capacity: 130 lb / 65 kg
- Graduation: 1 oz / 20 g
- Tray Size: 24-1/8" (w) x 14-7/8" (d) x 2-1/2" (h)
- Easy to read rotation poise bars for dual read in lb and kg - engraved numerals
- Easy to clean plastic tray
- Measuring Tape: 0" - 23" / 0 cm - 58 cm
- 2 Year Limited Warranty