




## Advantages

- ✓ **Sanitary Plastic Weighing Surface**  
The easy-to-clean plastic seat is ideal for sterile environments
- ✓ **Adjustable Unit Measurement Bar**  
The mechanical rotating poise bars allow for easy reading of weight in lb or kg
- ✓ **Engraved Numerals**  
Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
1524KL	1	20.75" x 45.75" x 18.5" / 5.5	44 lb	892076002183



**Specifications**

- Capacity: 130 lb / 65 kg
- Graduation: 1 oz / 20 g
- Seat Size: 14-1/2" (w) x 17-1/4" (d) x 18-1/4" (h)
- Easy to read rotation poise bars for dual read in lb and kg - engraved numerals
- Easy to clean plastic seat
- 2 Year Limited Warranty