

Health o meter®

P r o f e s s i o n a l

EMR



Advantages

- EMR/EHR Connectivity Featuring USB Port**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB
- Hospital Grade Construction**
Durable steel structure with non-slip platform enables patients to be weighed safely, accurately and comfortably
- “Live” Handrails and 180° Swivel Display**
Handrails may be used during weighing for patient stability and swivel display allows for patient privacy and ease of use
- Value-Added Features**
Operating functions: BMI, Tare, Zero Out, Hold, Live Weight

Health o meter®
P r o f e s s i o n a l
www.homscales.com

See reverse side for additional specifications

“LIVE” Handrails
for patient stability



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2101KL	1			892076002282



Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- 3/4” LCD display, with 180° swivel for ease of viewing
- Functions: LB / KG Conversion, LB / KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- Handrails may be used during weighing for patient stability
- Platform Size: 26” (w) x 22” (d)
- (2) Wheels for ease of movement
- 120V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty



2101KL: Reduced or minimum foam in packaging / recyclable cardboard packaging.