

*Low Acquisition Cost, Exclusive Product Technology
and Industry Leading Product Warranty*



Advantages

- Innovative Technology**
Motion sensing technology compensates for baby's movement and high-resolution readings provide more accurate weight measurements
- Industry Leading 3 Year Warranty**
Protects you against damage and repair
- Value-Added Features**
Operating functions (Tare, Zero Out, LB/KG Lock Out, Auto Zero, Auto Off) offered at no additional cost
- Flexible Functionality**
Plastic easy to clean tray ideal for sterile environments and heavy-duty built-in handles allow for easy portability

Healthometer®
P r o f e s s i o n a l

Healthometer®
P r o f e s s i o n a l
www.homscales.com

See reverse side for additional specifications

HIGH RESOLUTION
PORTABLE SCALE



NON-CONDUCTIVE
WHEELS

2210CART

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2210KL	1	25.75" x 19.5" x 7" / 2.03	30 lb	892076002329
2210CART	1	2 Cartons: 26.5" x 18.5" x 9" / 2.55	23 lb	
		22" x 17.75" x 2.5" / 0.56	17 lb	

2210KL



2210KL Specifications

- Capacity: 45 lb / 20 kg
- Resolution: 0.1 oz / (0 -11300 g / 1 g; >11300 g / 5 g)
- 1 gram accuracy standard
- 7/8" LCD display
- Functions: LB / GM Conversion, GM / LB Lock Out, Zero Out / Tare, Reweigh / Recall, Auto Zero, Auto Off, Motion Detection
- Removable curved clear acrylic tray
- Tray Size: 24-1/16" (w) x 15-5/16" (d) x 3-1/2" (h)
- Measuring Tape: 0" - 23" / 0 cm - 60 cm
- Heavy duty handles for ease of moving
- 6 C-cell batteries included
- 3 Year Limited Warranty



2210KL: No foam in packing and recyclable cardboard packaging

2210CART Specifications

- Top shelf has handle for easy grip
- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h)
- Lockable Cabinet: 23" (w) x 16-3/4" (d) x 17" (h)
- Storage Shelf is easy to clean, durable blue plastic
- Use with pediatric scales: 2210KL, 553KL, 549KL
- 1 Year Limited Warranty