


Remote Display can be mounted on a wall or placed on a table



Advantages

- Connectivity via Serial Port**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via Serial Port for fast, easy and accurate data transfer
- Durable Construction**
The scale has an easy-to-clean stainless steel solid construction and durable, non-skid mat and can withstand frequent use
- 550 lb Capacity Floor Scale**
The high capacity of this floor scale provides more versatile use
- AC Adapter Included**
Scale is powered by a 120V AC adaptor (included) or 4 AA batteries (not included) which allows for the convenience of dual power sources
- Motion Detection**
The advanced technology of this scale can accurately weigh patient while compensating for movement
- Ideal for Weighing Laundry**
The 42-1/8" x 21-3/4" platform conveniently accommodates standard medical linen carts

Remote Display can be mounted on a wall or placed on a table



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2842KL	1	49.5" x 24.25" x 4" / 2.78	65 lb	892076002480



Specifications

- Capacity: 550 lb / 250 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via Serial Port
- 7/8" LCD display
- Functions: LB / KG Conversion, Zero Out / Tare, Motion Detection
- Platform Size: 42-1/8" (w) x 21-3/4" (d)
- 120V adapter included, 4 AA batteries (not included)
- Can be used for laundry
- 3 Year Limited Warranty



Reduced or minimum foam in packaging / recyclable cardboard packaging.