



Advantages

- ✓ **Adjustable Arm & Foot Rest**
The flip-arm and foot rest allows for easier patient access and patient comfort
- ✓ **Sanitary Plastic Weighing Surface**
The easy-to-clean plastic seat is ideal for sterile environments.
- ✓ **Wheels**
The two back wheels provide easy mobility and allows for convenient placement



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
445KL	1	33.25" x 23.5" x 41.13" / 18.04	81 lb	892076002046



Specifications

- Capacity: 440 lb / 200 kg
- Graduation: 1/4 lb / 100 g
- Seat Size: 16-1/2" x 18-1/2" x 17-1/8"
- (2) Wheels for ease of movement, flip arm rest and foot rest
- Easy to clean plastic seat
- 2 Year Limited Warranty



Reduced or minimum foam in packaging / recyclable cardboard packaging.