



Remote Display can be mounted on a wall or placed on a table

Advantages

- EMR/EHR Connectivity via USB**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer
- Body Mass Index (BMI) Capability**
Scale functions include the ability to automatically calculate the patient's Body Mass Index (BMI)
- Large Platform Size**
Oversized 17-3/4" (w) x 14" (d) platform provides a safer and more accommodating weighing surface
- Value-Added Features**
Operating functions (BMI, Tare, Zero Out, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Remote Display**
The remote display provides convenience with the ability to be mounted on a wall or placed on a table



Remote Display can be mounted on a wall or placed on a table

PORTROD
Optional Wall Mounted
Height Rod Available



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
498KL	1	25.75" x 17.25" x 5.275" / 1.38	16 lb	892076002718



Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- 1" LCD display
- Remote Display can be mounted on a wall or placed on a table
- Functions: LB / KG conversion, LB /KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size: 17-3/4" (w) x 14" (d)
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- Optional Wall Mounted Height Rod available (item # PORTROD)
- 2 Year Limited Warranty



Carrying Case for 498KL
Item# 64771