



Advantages

- ✓ **EMR/EHR Connectivity via USB**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer
- ✓ **Value-Added Features**
Operating functions (LB/KG Lockout, Tare, Zero Out, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **1" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **Convenient Weighing Seat**
Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
524KL	1	17.13" x 27" x 20.5" / 5.49	26 lb	892076002534



Specifications

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb 0.5 oz 0 - 9 kg / 5g; 9 - 23 kg / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG conversion, LB /KG lock out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Seat Size: 18" (w) x 14-1/8" (d) x 18-1/2" (h)
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty