



Advantages

- ✓ **EMR/EHR Connectivity via USB**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer
- ✓ **Value-Added Features**
Operating functions (LB/KG Lockout, Tare, Zero Out, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Motion Detection**
The advanced technology of this scale can accurately weigh patient while compensating for movement
- ✓ **Sanitary Plastic Weighing Surface**
The easy-to-clean plastic tray is ideal for sterile environments
- ✓ **Built-In Measuring Tape**
Easy-to-read measuring tape allows user to measure both weight and height simultaneously



2210CART
Rolling Scale Cart for 553KL

- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h)
- Lockable cabinet: 23" (w) x 16-3/4" (d) x 17" (h)
- Storage Shelf
- 1 Year Limited Warranty

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
553KL	1	27.75" x 19.13" x 5.25" / 1.61	12 lb	892076002329

Specifications



- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Tray Size: 25" (w) x 13-7/8" (d) x 2" (h)
- Measuring tape: 0" - 24" / 0 cm - 58 cm
- 6 AA batteries included, optional power adapter (order # ADPT30)
- Optional carrying case (order # 553CASE)
- Optional Rolling Scale Cart (order # 2210CART)
- 2 Year Limited Warranty



Carrying Case for 553KL
Item # 553CASE