

- ✓ 100% Digital
- ✓ Stable, No Shaking



Advantages

- ✓ **Value-Added Features**
Operating functions (Tare, Zero Out, Auto Off, LB/KG Lockout) offered at no additional cost
- ✓ **Convenient Power Source**
Scale is powered by a 6V rechargeable battery and includes battery charger
- ✓ **1-1/4" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **100% Digital Weighing System**
This 100% digital scale is easy-to-use and efficiently provides accurate weight measurements
- ✓ **Adjustable Arm & Foot Rest**
The flip-arm and foot rest allows for easier patient access and patient comfort
- ✓ **Sanitary Plastic Weighing Surface**
The easy-to-clean plastic seat is ideal for sterile environments
- ✓ **Sturdy Construction**
The sturdy construction of this scale provides a stable base that does not shake

Health o meter®
P r o f e s s i o n a l



- 100% Digital
- Stable, No Shaking



594KL Display

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
594KL	1	32.25" x 23.5" x 41.13" / 18.04	81 lb	892076002688



Specifications

- Capacity: 600 lb / 280 kg
- Resolution: 0.2 lb / 0.1 kg
- 1-1/4" LCD display
- Functions: LB / KG Conversion, LB / KG lock out, Zero Out / Tare, Auto Zero, Auto Off
- (2) Wheels for ease of movement, flip arm rest and foot rest
- Seat Size: 18-1/4" x 14-1/2" x 17-1/2"
- Easy to clean plastic seat
- Charger included, 6 Volt rechargeable battery included
- 2 Year Limited Warranty



No foam in packing and recyclable cardboard packaging