



8320KL Open

Advantages

- Sanitary Plastic Weighing Surface**
The easy-to-clean plastic tray is ideal for sterile environments
- Motion Detection**
The advanced technology of this scale can accurately weigh patient while compensating for movement
- Tare Function**
This convenient function allows the user to automatically subtract the weight of blanket or diaper from infant's weight measurement
- Value-Added Features**
Operating functions (Tare, Zero Out, Auto Zero, Auto Off) included at no additional cost
- Easy Portability**
This scale is ideal for traveling nurses or small healthcare facilities with its sturdy construction, built-in carrying handles, easy to clean surface, optional carrying case and ability to fold in half



8320KL Open

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
8320KL	1	22" x 18" x 14.8" / 3.39	24.69 LB	892076002817

Specifications

- Capacity: 33 lb / 15 kg
- Resolution: 0.5 oz / 10 g
- 3/4" LCD display
- Functions: LB / KG Conversion, Zero Out / Tare, Auto Zero, Auto Off
- Tray Size: 23-1/2" (w) x 11" (d) x 1" (h)
- 2 AAA batteries included
- Optional carrying case (order # 64771)
- 2 Year Limited Warranty



Carrying Case for 8320KL
Item # 64771