

Powered by

**VALHALLA**  
SCIENTIFIC

Body Composition Scale

## ADVANTAGES

### 15 Second Measurement Time

- ✓ On screen display delivers rapid access to weight, body fat %, body water %, lean mass, BMR & BMI with a full illustrated printout in 60 seconds.

### Large Platform

- ✓ The spacious platform allows patients to find a comfortable and secure stance.

### Few Pretest Requirements

- ✓ Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

### Safe to Use

- ✓ This scale is FDA registered.

### High Capacity

- ✓ The scale features an industry leading 1000 lb capacity, that provides the flexibility and accuracy to weigh any patient.

### Patient Report Templates\*

- ✓ The G6 Series scales offer users 3 graphic printout formats to use with their patients: Adult/Bariatric, Pediatric and Fitness.

\* User must connect scale to a printer to provide personalized graphic printout. Printer is optional and sold separately, not available through Health o meter® Professional scales.

Powered by

**VALHALLA**  
SCIENTIFIC



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
BCS-G6	46" x 19" x 10" / 10	37 lb	854896003702

**Specifications**

- Capacity: 1000 lb / 453.5 Kg
- Resolution: 0.1 lb / 0.1 kg
- Platform Size: 14" (w) x 16" (d) x 1.875" (h)
- Display: 6-window display (.5" digits) includes weight, % body fat, % body water, % muscle mass, BMI, REE
- Power Source: 110-220V power adapter (included)
- Printer Connectivity: USB or WiFi (optional printer, sold separately, not available through Health o meter® Professional scales)
- 1 Year Limited Warranty

**Bio-Impedance Specifications**

- Electrode Method: Tetrapolar System
- Frequency: 50kHz
- Test Current: 500µA
- Conversion Rate: 750msec

Model	Description	UPC Code
BCS-G6-adult	Body Comp Scale - Adult Body	854896003702
BCS-G6-limbs	Body Comp Scale - Adult Including Limbs	854896003719
BCS-G6-DUO	Body Comp Scale - Adult and Pediatric	854896003733
BCS-G61-upper	Body Comp Scale - Upper Body Only	854896003726
IPO-ADU	Patient Report Template - Adult (5-page report, 25 reports/package)	854896003740
IPO-PED	Patient Report Template - Pediatric (7-page report, 25 reports/package)	854896003757
IPO-ATH	Patient Report Template - Athletic (3-page report, 50 reports/package)	854896003764
BCS-WIRELESS	Body Comp Scale Wireless Adapter	854896003771

Body Composition Scale

