



Harmar®

AL301XL / AL301XL HD Extra Large Fusion Lifts

Extra Large Fusion Lifts

The AL301XL's larger deck handles wider or longer scooters and power chairs. It's the ideal choice to gain 8.5" extra width or meet the need for lengths longer than 46." Available in 350 lbs. or heavy-duty 400 lbs. capacity models.

Features & Benefits

- Easy, one-switch operation
- Self-tensioning Q'Straint retractors and front and rear wheel cradles; no chair modifications required
- 36" x 39" aluminum deck; adjustable cradle extends up to 53"
- Manual crank backup and license plate mount included
- Swing-Away option available



HARMAR EXCLUSIVE

Specifications	AL301XL	AL301XL HD Heavy Duty
Accommodates	Larger power chairs & scooters	Larger power chairs & scooters
Platform	36" x 39" aluminum deck; adjustable cradle extends up to 53"	36" x 39" aluminum deck; adjustable cradle extends up to 53"
Lifting Capacity	350 lbs. (160 kg)	400 lbs (181 kg)
Installed Weight	100 lbs. (45 kg)	102 lbs. (46 kg)
Securement	Self-tensioning Q'Straint retractors and front and rear wheel cradle	Self-tensioning Q'Straint retractors and front and rear wheel cradle
Hitch	Class III or greater	Class III or greater
Options	Battery pack, Swing-Away, Power chair or scooter covers	Battery pack, Swing-Away, Power chair or scooter covers
Warranty	3-year transferable	3-year transferable



Harmar® 800-833-0478 | www.harmar.com
America's Lift Leader™



© 2014 Harmar MKT-Harmar SS-AL301XL 12Feb14 Part Number: MKT-000000 Rev A